


How to Process GCMS Data

After you run the GCMS to collect data for your sample, you need to get the mass spectrum of your product.

You will use a copy of the GCMS software that *is not attached to the instrument* to print your data.



Load the CASPiE Instrument Site



- To open the GCMS data processing software, launch your web browser and navigate to the CASPiE instrument website at <https://instruments.caspie.org>.
- Enter your account information to login.

The image shows a login interface for the CASPiE Instrument Access system. At the top, the CASPiE logo is on the left and the text "Instrument Access" is on the right. Below this is a grey login box with the heading "Welcome CASPiE Users" and the instruction "Login to access your CASPiE resources." Inside the box are two input fields: "User name:" with the text "jmsmith" and "Password:" with a masked password of eight dots. A "Log On" button is located at the bottom right of the login box. Below the login box, a paragraph of text provides instructions and contact information for users having trouble logging in.

CASPiE Instrument Access

Welcome CASPiE Users

Login to access your CASPiE resources.

User name: jmsmith

Password:

Log On

Please enter your account credentials then click Log On to get access to your reserved instrument time. If you have trouble logging in, contact Debora Steffen by [email](#) or phone (765-494-4959).

The GCMS Data Processing Software




You will use a copy of the GCMS software that *is not attached to the instrument* to process your data files.

The stand alone version of the software called '**GCMS Data Processing**'.

You will always have access to the GCMS data processing software during a CASPiE module. You will only have access to the GCMS instrument for remote control during your scheduled hours.

Important Note:
You must have the Citrix client installed on your computer to proceed!




Launch GCMS Data Processing Software



- You should see the **GCMS Data Processing** icon in the Applications box. Launch the software by clicking once on the icon (a). If the **GCMS Data Processing** icon is not present ask your TA for help.
- A login window will appear (b). Use the user ID 'Student.' *No password is required.*


(a)



Open Data File


- The data analysis application that opens is called **GCMS Postrun Analysis**.
- In the Data Explorer window, double-click your data file to open the file.

Data Explorer window




Data Display

- The data display has three windows, two chromatogram windows and the mass spectrum window.



Chromatogram Windows


- The bottom chromatogram window can display an expanded view of the chromatogram.
- In the top window click and drag around the region of interest. The expanded view of the selected region will be shown in the lower window.
- You can adjust the y-scale of the chromatogram using the tools shown.
- Double click on any point in the chromatogram to display the mass spectrum at that retention time.




Integrate Peaks

- In order to produce a quality spectrum of your peak(s), you must integrate the chromatographic peaks first.
- Select 'Qualitative| Peak Integrate| TIC' (a) from the menu.
- The Qualitative Parameters window will open (b).
- Change the minimum area to 100,000 and click OK. Keep the other default values.

(a)





(b)



Optimize Integration


- If there are too many peaks marked (lots of red arrows) change the minimum integration area to a larger number.
- Select 'Qualitative| Peak Integrate| TIC' again and change the minimum area by one order of magnitude and re-integrate by clicking OK.
- Repeat this process, increasing the minimum peak area until you get only the peaks of interest marked by red arrows. **Make sure all the large peaks are marked.**
- How many peaks you mark is up to you. Generally, the very small peaks are not of interest (insufficient concentration to be significant).




Launch Report Template

- To print the chromatogram and peak spectra you must create a report.
- Click on the 'Report' icon on the left (a).
- The report template file will be displayed in the Data Explorer window (b).
- Double click on the 'Antiviral' file to open the report template (b).

(a)




(b)




Print Report

- To print the report select File |Print (a).
- Your local printer and an Adobe PDF printer should be listed in the print window (b).
- You can print directly to your printer or create a PDF file.
- If you want to create a PDF file, select the Adobe PDF printer and follow the directions on the page 14.

(a)




(b)




The Report

- The report includes sample information, the chromatogram, and mass spectra of all integrated peaks.



The Report



- The printed mass spectrum is the peak average minus a background spectrum.



Optional PDF Printing


- If you want to create a PDF file, select File| Print and select the 'Adobe PDF' printer. Then click on the PDF printer properties button (a).
- On the 'Adobe PDF Settings' tab click on the 'Browse' button next to the 'Adobe PDF Output Folder' box. Select the Z:\ drive folder. This is your folder on the CASPiE file server.
- The rest of the settings should be as shown.

(b)



Close Windows


When you have finished printing exit the GCMS software.



External Access to Data Files

- You can access your data files without the GCMS software. To do this use your browser to navigate to <https://data.caspie.org>.
- Login with your user name and password (a) to display the data in your CASPIE folder.
- You can download any pdf file you may have created by right clicking on the file and selecting 'Save Target As....' (b).

(a)



(b)

